



cap-ONE
mask

Ensure Quality of Care under Sedation

- Accurate respiratory monitoring with
optimal oxygen supply -

Efficient Monitoring for Safer Pediatric Sedation

— Don't miss signs of respiratory depression —

Pediatric sedation may involve serious risks including hypoventilation, apnea, airway obstruction and cardiopulmonary dysfunction¹). Also, it has been shown that depth of sedation can be deeper than the intended level in pediatric patients²). For safer patient management, it is essential to closely observe patients while administering sedation. According to an analysis of adverse events related to pediatric sedation in the US, the most commonly observed adverse event is respiratory depression. It also reports that these adverse events can lead to serious incidents including death and permanent dysfunction in the absence of appropriate monitoring³).

— Fast, accurate detection of respiratory depression while oxygen is supplied —

Nihon Kohden has developed an extremely small mainstream CO₂ sensor (cap-ONE) that can accurately measure oral as well as nasal exhaled CO₂ in non-intubated patients⁴). After further development and evolution, the cap-ONE is now available as the world's smallest mainstream CO₂ sensor weighing only 4 g with an original open face oxygen mask (cap-ONE mask). Cap-ONE mask is easily used for pediatric patients who are receiving supplemental oxygen, and the cap-ONE CO₂ sensor reliably detects respiratory depression including apnea and hypoventilation while oxygen is supplied (Fig. 1).



cap-ONE
ORAL NASAL EXPIRATION

Using this ultra-small mainstream CO₂ sensor together with cap-ONE mask, the sensor can accurately catch oral and nasal exhaled CO₂ and reliably detect respiratory variation.




Figure 1.

The exhaled CO₂ waveform (capnogram) simulated by a pediatric model mannequin and spontaneous breathing lung simulator in the presence of oxygen supplement. The cap-ONE sensor reliably detects exhaled CO₂ during oxygen administration and presents a quality capnogram which agrees with a reference capnogram. (From Reference 5, with permission)

— Gentle oxygen supply from a mask with a wide opening —



cap-ONE
mask

Gentle stream of oxygen from both sides of the mask

Consistent and stable oxygen concentration

Less rebreathing of exhaled CO₂ even with low oxygen flow

Lighter, freer, less oppressive feeling

It is recommended that oxygen flow should be more than 6L/min when using an oxygen mask in order to prevent the patient from rebreathing exhaled gas which is retained in the mask⁶⁾. Cap-ONE mask has a wide opening in the front to let oral and nasal exhaled gas out of the mask. This design prevents exhaled gas from remaining in the mask even when oxygen flow is low (Fig. 2). Also, well dispersed oxygen is supplied from the sides of cap-ONE mask, so oxygen concentration inside the mask remains consistent and stable⁷⁾.




Figure 2.

Comparison of CO₂ partial pressure in exhaled gas while administering oxygen using a conventional oxygen mask and cap-ONE mask (pediatric patients aged 1 to 8, CO₂ partial pressure measured 5 minutes after the oxygen administration). Point A shows retained exhaled gas and increased CO₂ partial pressure in the conventional mask. Point B shows that CO₂ partial pressure remains low and rebreathing of exhaled gas is reduced. (From Reference 7, with permission)

Limitations of the Pulse Oximeter

— Desaturation can be delayed during oxygen administration —

The pulse oximeter is widely used to monitor vital signs in various clinical settings since it provides oxygen saturation and pulse rate non-invasively and continuously with just a probe on the patient’s finger. However, it is important to understand the limitations of pulse oximeters especially when using it as a respiratory monitoring tool. It has been shown that when fraction of inspiratory oxygen (FiO_2) is high, saturation level remains high and does not decrease promptly even when pulmonary function is deteriorated⁸⁾ (Fig. 3). This can lead to failure to detect respiratory depression and hypoventilation for a long time. There has been a reported case where a post-operative patient who received high concentration oxygen developed serious and fatal hypercapnia and acidosis in the recovery room⁹⁾. Monitoring respiration only with a pulse oximeter can have potential risks. It is important to understand that the pulse oximeter can be an excellent index of oxygenation only when a patient breathes room air¹⁰⁾.

In June 2011, the Anesthesia Patient Safety Foundation held a meeting titled “Essential Monitoring Strategies to Detect Clinically Significant Drug-Induced Respiratory”. Different monitoring methods were compared (Table 1) and capnography (monitoring of exhaled CO_2) was reported to be the most suitable method for detecting respiratory depression during oxygen administration¹⁰⁾.




Figure 3. Change of oxygen saturation plotted as a function of increasing intrapulmonary shunt from 0% at a rate of 2%/min. With fraction of inspiratory oxygen (FiO_2) of 1.0, 0.6, and 0.3, it takes 16, 11, and 6 minutes for oxygen saturation (SaO_2) to decrease from 100% to about 97%, respectively. (From Reference 8 with permission, partly modified)

Monitoring method	Sensitivity	Specificity	Reliability	Response time
$P_{et}CO_2$ (intubated)	High	High	High	Fast
$P_{et}CO_2$ (not intubated)	High	Moderate-High	Moderate	Fast
S_pO_2 (no O_2 supplement)	High	Moderate-High	High	Fast
S_pO_2 (with O_2 supplement)	Moderate	Moderate	High	Slow
Clinical assessment by skilled clinician	Moderate	Moderate-High	Moderate	Slow
Thoracic impedance	Low-Moderate	Low	Low	Moderate

Table 1. Monitoring methods for detecting opioid induced respiratory depression in the postoperative period. (From Reference 10, with permission)

Exhaled CO₂ Monitoring Recommended by Anesthesia Guidelines

Reflecting serious concerns about pediatric sedation outside the operating room, the American Academy of Pediatrics and the American Academy of Pediatric Dentistry updated their “Guidelines for Monitoring and Management of Pediatric Patients During and After Sedation for Diagnostic and Therapeutic Procedures”¹⁾ in 2006. This guideline describes recommendations regarding monitoring of exhaled CO₂ as well as monitoring of common vital signs including heart rate and oxygen saturation. Also, the “Standards for Basic Anesthetic Monitoring” (revised July 2010 and effective July 2011) mandated that the adequacy of ventilation should be evaluated by exhaled CO₂ in moderate to severe sedation¹¹⁾. Now, capnography is strongly recommended as a means to continuously assess patients’ respiratory condition during procedures using sedation and anesthesia which may be associated with risks of respiratory depression.

cap-ONE Product Line



cap-ONE CO₂ sensor
P909 (TG-970P)

For non-intubated pediatric patients



cap-ONE mask
V933 (YG-232T)

For child weighing
20-40kg



cap-ONE mask
V935 (YG-242T)

For child weighing
7-20kg

For intubated patients



Airway adapter
R805 (YG-211T) **A**

For adult and child weighing over 7kg

R806 (YG-213T) **B**

For infant weighing 2-7kg

cap-ONE Information Vol.2



References

- 1) Cote CJ et al. Guidelines for Monitoring and Management of Pediatric Patients During and After Sedation for Diagnostic and Therapeutic Procedures: An Update (American Academy of Pediatrics). PEDIATRICS 2006; Vol.118, No.6: 2587-2602.
- 2) Motas D et al. Depth of consciousness and deep sedation attained in children as administered by non-anaesthesiologists in a children's hospital. Paediatr Anaesth 2004; Vol.14: 256-260.
- 3) Cote CJ et al. Adverse Sedation Events in Pediatrics: A Critical Incident Analysis of Contributing Factors. PEDIATRICS 2000; Vol.105, No.4: 805-814.
- 4) Kasuya Y et al. Accuracy of Postoperative End-tidal PCO₂ Measurements with Mainstream and Sidestream Capnography in Non-obese Patients and in Obese Patients with and without Obstructive Sleep Apnea. Anesthesiology 2009; Vol.11, No.3: 609-615.
- 5) Takatori F et al. A Novel Mainstream Capnometer System for Non-intubated Pediatric Patients Requiring Oxygen Administration. IEEE Eng Med Biol Soc. 2011; 2011; 1189-1192.
- 6) American Heart Association. Respiratory Management Resources 2006.
- 7) Abe. S et al. 2010. Japanese Society of Pediatric Intensive and Critical Care. A novel method for mainstream capnometer. Presentation. November.
- 8) Downs JB, Smith RA. Increased inspired oxygen concentration may delay diagnosis and treatment of significant deterioration in pulmonary function. Crit Care Med.1999 Dec;27(12):2844-6.
- 9) Ayas N et al. Unrecognized Severe Postoperative Hypercapnia: A Case of Apneic Oxygenation. Mayo Clinic Proceedings; Vol.73, No.1: 51-54.
- 10) Weinger MB et al. No Patient Shall Be Harmed By Opioid-Induced Respiratory Depression. APSF Newsletter 2011; Vol.26, No.2: 21-28.
- 11) American Society of Anesthesiologists. Standards for Basic Anesthetic Monitoring. <http://www.asahq.org/For-Members/Standards-Guidelines-and-Statements.aspx>



NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku, Tokyo 161-8560, Japan
Phone +81 (3) 5996-8036 Fax +81 (3) 5996-8100
www.nihonkohden.com