**ITALIANO**

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Nome: Isotonac-3, Isotonac-4, Hemolynac-3, Hemolynac-3N, Hemolynac-5, Hemolynac-310, Hemolynac-510, Cleanac-3, Cleanac, Cleanac-710

**DESTINAZIONE D'USO**

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**BIBLIOGRAFIA**

1. Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard – Second Edition, CLSI Document H26-A2, Vol. 30, No. 14, 2010.
2. Health Care Financing Administration. Clinical Laboratory Improvement Amendments of 1998; Final Rule. Federal Register, Vol. 57, [42 CFR § 493.1217, 493.1218], 1992.

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**NIHON KOHDEN FIRENZE S.r.l**

Via Torta, 72/74 - 50019 Sesto Fiorentino-Firenze, Italy  
Tel +39 055 3045 1 - Fax +39 055 308548

**ENGLISH**

**PRODUCT IDENTIFICATION**

name: Isotonac-3, Isotonac-4, Hemolynac-3, Hemolynac-3N, Hemolynac-5, Hemolynac-310, Hemolynac-510, Cleanac-3, Cleanac, Cleanac-710

**INTENDED USE**

For in vitro diagnostic use only. Reagents to be used on NIHON KOHDEN analyzers in combination as indicated in Table 1.

Isotonac-4, Hemolynac-310, Hemolynac-510 and Cleanac-710 are only intended for human blood. In regard to the intended use, please refer to Table 2.

**STORAGE AND SHELF LIFE AFTER FIRST OPENING**

Store between 1 and 30 °C. Check reagent integrity and expiry date, before installing the device. Once opened, the reagent must be used within the expiry date indicated in Table 4. When not in use, use the original cap to close the container except Hemolynac-310, Hemolynac-510 and Cleanac-710. Do not force.

**COMPOSITION**

Active ingredients

Ref. Table 3

**WARNING AND PRECAUTIONS**

Labelling information

Ref. Table 5.

Disposal information

Unused product and emptied containers should be disposed in accordance with local environmental regulations. Instrument wastes, as potential biohazardous, must be handled with care and disposed in accordance with local environmental regulations.

**SPECIMEN**

The devices have been designed to be used with whole blood samples. Refer to Operator's Manual for sample collection and storage conditions.

**PROCEDURE**

Refer to Operator's Manual for reagents installation and priming, and for instrument operation. When installing a new batch of diluent or lysing reagent, check the background according to Operator's Manual. Refer to Operator's Manual for working temperature.

Use each reagent in its own packaging (do not transfer). Follow GLP recommendations. If the instrument connections were not compatible, work according to Good Laboratory Practice (GLP), in order to avoid any kind of contamination, either from external environment to the device, or from the device to the external environment.

**METHODOLOGY AND CALCULATION OF ANALYTICAL RESULTS**

Refer to Operator's Manual for information on methodology and calculation of analytical results.

**INTERNAL QUALITY CONTROL**

Background verification

Before running testing samples check background count according to Operator's Manual, and verify that all the values for each parameter are within Operator's Manual limits. If not, follow Operator's Manual prescriptions. In case of problems attributed to one of the reagents, install a new unit of reagent and verify the background again. Do not run testing samples until the background is within the specified limits.

Calibration and calibration verification procedures

Instrument calibration should be verified in the following conditions:

1. after first installation
  2. at least every six months (or according to manufacturer's requirements for calibration frequency, if more stringent)
  3. after a major repair or a major change (replacement of critical parts that may influence test performance, software upgrade, complete change of reagents)
  4. when quality control data reflect a significant change of baseline accuracy (drift), or an out-of-control condition
  5. in case of replacement of the lysing reagent Hemolynac-3 with the lysing reagent Hemolynac-3N and vice versa
- Calibration procedures should be preceded by a preventive maintenance (verification of baseline accuracy and precision: ref. CLSI document H26-A2). Calibration and calibration verification must be performed using appropriate calibration material (fresh whole blood samples which have been analyzed by reference methods, or appropriate calibrator), following Operator's Manual prescriptions. When calibration is required, verify instrument calibration immediately after calibrating, according to Operator's Manual instructions (refer also to CLSI document H26-A2).

Daily Quality Control Procedures

In order to assure the accuracy and reliability of test results, it is recommended to perform quality control procedures, according to Operator's Manual instructions, using appropriate control material (retained patient specimens or stabilized control blood), immediately before, during, or immediately after each analysis run (intended as an analysis sequence). It is recommended to run two different levels of stabilized control blood (refer to CLSI document H26-A2).

**BIBLIOGRAPHY**

1. Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard – Second Edition, CLSI Document H26-A2, Vol. 30, No. 14, 2010.
2. Health Care Financing Administration. Clinical Laboratory Improvement Amendments of 1998; Final Rule. Federal Register, Vol. 57, [42 CFR § 493.1217, 493.1218], 1992.

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