



# SSD Expansion Board

## YS-116P0

### General

The YS-116P0 SSD expansion board is for the CSM-1502, CSM-1701 or CSM-1702 bedside monitor (CU-152R, CU-171R or CU-172R core unit). By installing the SSD expansion board in the monitor, a large capacity of storage can be obtained so that more full disclosure waveforms can be saved or more applications can be installed in the monitor.

### Safety Information

For safety information, refer to the Operator's Manual or Installation Guide of the device using this product.

This Safety and Performance Information is an extract from the general and safety information sections of the most recent edition of Operator's Manual or Installation Guide. Therefore, the contents of your Operator's Manual or Installation Guide may differ from those of this Safety and Performance Information. For detailed operating procedures, follow the instructions of your Operator's Manual or Installation Guide.



Manufacturer

**NIHON KOHDEN CORPORATION**  
1-31-4 Nishiochiai, Shinjuku-ku,  
Tokyo 161-8560, Japan  
Phone +81 3-5996-8041  
<https://www.nihonkohden.com/>



European Representative

**NIHON KOHDEN EUROPE GmbH**  
Raiffeisenstrasse 10, D-61191 Rosbach, Germany  
Phone +49 6003-827-0 Fax +49 6003-827-599