

MEK-5DN HEMATOLOGY CONTROL
for WBC 5 part differential 3.0mL

ASSAY SHEET

LOT	MEK5D242N
	2024-05-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±0.8	8.1	8.1	7.9	8.11
NE%	%	±15.0	44.7	44.7	44.7	44.46
LY%	%	±15.0	37.0	37.0	37.0	39.83
MO%	%	±5.0	13.0	13.0	13.0	10.68
EO%	%	±5.0	3.1	3.1	3.1	3.57
BA%	%	±4.0	2.2	2.2	2.2	1.46
NE	10 ³ /μL	±1.5	3.5	3.5	3.6	3.60
LY	10 ³ /μL	±1.5	3.0	3.0	2.9	3.23
MO	10 ³ /μL	±0.5	1.1	1.1	1.0	0.87
EO	10 ³ /μL	±0.5	0.3	0.3	0.2	0.29
BA	10 ³ /μL	±0.4	0.2	0.2	0.2	0.12
RBC	10 ⁶ /μL	±0.20	4.53	4.45	4.51	4.60
HGB	g/dL	±0.5	13.6	13.5	13.4	13.77
HCT	%	±3.0	42.4	40.5	40.0	47.6
MCV	fL	±6.0	93.7	91.1	88.6	103.4
MCH	pg	±2.8	30.0	30.3	29.7	29.9
MCHC	g/dL	±3.0	32.1	33.3	33.5	28.9
RDW-CV *	%	±4.0	13.8	14.6	14.0	14.2
PLT	10 ³ /μL	±50	276	293	269	271.7
MPV	fL	±2.0	3.8	5.3	4.0	7.0

Manual Registration	
Lot	242
Level	Normal
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901691N