


MEK-3DH HEMATOLOGY CONTROL 2.0mL


ASSAY SHEET

LOT	B245H
	2024-08-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order. If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±2.2	21.0	21.0	20.57
RBC	10 ⁶ /μL	±0.25	5.95	5.95	5.95
HGB	g/dL	±0.7	19.6	19.5	19.54
HCT	%	±4.0	54.4	54.4	62.1
MCV	fL	±7.0	91.4	91.4	104.4
MCH	pg	±2.8	32.9	32.8	32.8
MCHC	g/dL	±3.0	36.0	35.8	31.5
PLT	10 ³ /μL	±70	543	543	546.9
MPV	fL	±3.0	5.6	8.1	8.5
RDW-CV *	%	±3.0	16.1	16.1	17.6
LY%	%	±10.0	17.3	17.3	13.28
MO%	%	±5.0	2.1	2.1	2.11
GR%	%	±10.0	80.6	80.6	84.61
LY	10 ³ /μL	±1.8	3.6	3.6	2.73
MO	10 ³ /μL	±0.7	0.4	0.4	0.43
GR	10 ³ /μL	±2.0	17.0	17.0	17.41

Manual Registration	
Lot	245
Level	High
Control	MEK-3D

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



MEK-3DH



B245H



2024-08-05



0694-901686P