


MEK-3DL HEMATOLOGY CONTROL 2.0mL

ASSAY SHEET

LOT	B245L
	2024-08-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order. If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.6	2.1	2.1	2.05
RBC	10 ⁶ /μL	±0.15	2.39	2.39	2.32
HGB	g/dL	±0.4	6.0	6.2	6.09
HCT	%	±2.0	17.8	17.8	19.7
MCV	fL	±9.0	74.5	74.5	84.9
MCH	pg	±2.4	25.1	25.9	26.3
MCHC	g/dL	±3.0	33.7	34.8	30.9
PLT	10 ³ /μL	±30	70	70	78.2
MPV	fL	±3.0	6.1	8.4	8.4
RDW-CV *	%	±3.0	17.5	17.5	15.8
LY%	%	±10.0	63.4	63.4	66.96
MO%	%	±5.0	8.4	8.4	4.51
GR%	%	±10.0	28.2	28.2	28.53
LY	10 ³ /μL	±0.5	1.3	1.3	1.37
MO	10 ³ /μL	±0.2	0.2	0.2	0.09
GR	10 ³ /μL	±0.5	0.6	0.6	0.59

MEK-1305
QC ESR
48±17 mm

Manual Registration	
Lot	245
Level	Low
Control	MEK-3D

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901684P