


**MEK-5DL** HEMATOLOGY CONTROL  
for WBC 5 part differential 3.0mL

**ASSAY SHEET**

LOT	MEK5D245L
	2024-08-05
OPEN	14DAYS

MEK-9100,9200  
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 <sup>3</sup> /μL	±0.6	3.1	3.1	3.2	3.14
NE%	%	±15.0	27.3	27.3	27.3	27.77
LY%	%	±15.0	55.0	55.0	55.0	57.46
MO%	%	±5.0	11.5	11.5	11.5	9.45
EO%	%	±5.0	3.5	3.5	3.5	3.94
BA%	%	±4.0	2.7	2.7	2.7	1.38
NE	10 <sup>3</sup> /μL	±0.5	0.8	0.8	0.8	0.88
LY	10 <sup>3</sup> /μL	±0.5	1.7	1.7	1.8	1.80
MO	10 <sup>3</sup> /μL	±0.2	0.4	0.4	0.4	0.30
EO	10 <sup>3</sup> /μL	±0.2	0.1	0.1	0.1	0.12
BA	10 <sup>3</sup> /μL	±0.2	0.1	0.1	0.1	0.04
RBC	10 <sup>6</sup> /μL	±0.15	2.66	2.64	2.67	2.61
HGB	g/dL	±0.4	6.7	6.7	6.6	6.73
HCT	%	±2.0	22.1	21.2	20.9	23.3
MCV	fL	±7.0	82.9	80.3	78.3	89.3
MCH	pg	±2.4	25.2	25.4	24.7	25.8
MCHC	g/dL	±3.0	30.3	31.6	31.6	28.9
RDW-CV *	%	±4.0	15.3	16.9	15.7	15.9
PLT	10 <sup>3</sup> /μL	±30	45	51	46	43.9
MPV	fL	±3.0	4.2	5.2	4.4	7.2

Manual Registration	
Lot	245
Level	Low
Control	MEK5D

\* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-9016900