

**MEK-5DN** HEMATOLOGY CONTROL  
for WBC 5 part differential 3.0mL

**ASSAY SHEET**

LOT	MEK5D245N
	2024-08-05
OPEN	14DAYS

MEK-9100,9200  
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 <sup>3</sup> /μL	±0.8	8.2	8.1	8.0	8.16
NE%	%	±15.0	48.8	48.8	48.8	47.73
LY%	%	±15.0	31.2	31.2	31.2	35.25
MO%	%	±5.0	13.6	13.6	13.6	10.87
EO%	%	±5.0	3.7	3.7	3.7	4.12
BA%	%	±4.0	2.7	2.7	2.7	2.03
NE	10 <sup>3</sup> /μL	±1.5	4.0	4.0	3.9	3.88
LY	10 <sup>3</sup> /μL	±1.5	2.6	2.5	2.5	2.88
MO	10 <sup>3</sup> /μL	±0.5	1.1	1.1	1.1	0.89
EO	10 <sup>3</sup> /μL	±0.5	0.3	0.3	0.3	0.34
BA	10 <sup>3</sup> /μL	±0.4	0.2	0.2	0.2	0.17
RBC	10 <sup>6</sup> /μL	±0.20	4.79	4.71	4.77	4.86
HGB	g/dL	±0.5	13.5	13.5	13.3	13.68
HCT	%	±3.0	42.5	40.7	40.1	47.5
MCV	fL	±6.0	88.8	86.4	84.0	97.7
MCH	pg	±2.8	28.2	28.7	27.9	28.1
MCHC	g/dL	±3.0	31.8	33.2	33.2	28.8
RDW-CV *	%	±4.0	14.8	16.1	15.1	16.0
PLT	10 <sup>3</sup> /μL	±50	246	261	240	224.7
MPV	fL	±2.0	3.2	5.1	3.5	6.7

Manual Registration	
Lot	245
Level	Normal
Control	MEK5D

\* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-9016910