

MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS245
	2024-06-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	9.1	9.1	9.3	9.2	9.0	9.35	9.17
RBC	10 ⁶ /μL	4.67	4.67	4.63	4.59	4.62	4.69	4.66
HGB	g/dL	13.8	14.0	13.7	13.6	13.5	13.90	13.87
HCT	%	40.3	40.3	41.5	40.3	39.1	46.1	45.7
MCV	fL	86.2	86.2	89.6	87.7	84.7	98.2	98.1
MCH	pg	29.6	30.0	29.6	29.6	29.2	29.6	29.8
MCHC	g/dL	34.2	34.7	33.0	33.7	34.5	30.2	30.4
PLT	10 ³ /μL	257	257	251	254	253	271.0	279.1
MPV	fL	6.8	8.8	6.8	7.6	6.8	9.3	9.2
RDW-CV *	%	16.4	16.4	15.1	16.5	15.4	16.2	16.2

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	167
FL	87
SD	103
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

