


<b>MEK-3DL</b>	HEMATOLOGY CONTROL	2.0mL
----------------	--------------------	-------

**ASSAY SHEET**

LOT	B248L
	2024-11-05
OPEN	14DAYS

<b>MEK-1301,1302,1305 Barcode Entry</b>
---

The barcodes may be read in any order.  
 If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
<b>WBC</b>	10 <sup>3</sup> /μL	±0.6	2.0	2.0	2.02
<b>RBC</b>	10 <sup>6</sup> /μL	±0.15	2.32	2.32	2.26
<b>HGB</b>	g/dL	±0.4	5.9	6.0	6.03
<b>HCT</b>	%	±2.0	18.1	18.1	20.0
<b>MCV</b>	fL	±9.0	77.9	77.9	88.5
<b>MCH</b>	pg	±2.4	25.4	25.9	26.7
<b>MCHC</b>	g/dL	±3.0	32.6	33.1	30.2
<b>PLT</b>	10 <sup>3</sup> /μL	±30	68	68	75.4
<b>MPV</b>	fL	±3.0	6.6	8.7	8.6
<b>RDW-CV *</b>	%	±3.0	17.2	17.2	15.9
<b>LY%</b>	%	±10.0	64.4	64.4	64.99
<b>MO%</b>	%	±5.0	10.1	10.1	6.37
<b>GR%</b>	%	±10.0	25.5	25.5	28.64
<b>LY</b>	10 <sup>3</sup> /μL	±0.5	1.3	1.3	1.31
<b>MO</b>	10 <sup>3</sup> /μL	±0.2	0.2	0.2	0.13
<b>GR</b>	10 <sup>3</sup> /μL	±0.5	0.5	0.5	0.58

<b>MEK-1305</b>
QC ESR
49±17 mm

Manual Registration	
Lot	248
Level	Low
Control	MEK-3D

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



**NIHON KOHDEN CORPORATION**  
 1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN



0694-901684Q