


MEK-5DL HEMATOLOGY CONTROL
for WBC 5 part differential 3.0mL

ASSAY SHEET

LOT	MEK5D248L
	2024-11-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±0.6	3.0	3.0	3.0	2.99
NE%	%	±15.0	27.9	27.9	27.9	26.75
LY%	%	±15.0	51.1	51.1	51.1	55.97
MO%	%	±5.0	13.3	13.3	13.3	10.42
EO%	%	±5.0	4.4	4.4	4.4	5.17
BA%	%	±4.0	3.3	3.3	3.3	1.69
NE	10 ³ /μL	±0.5	0.9	0.9	0.9	0.81
LY	10 ³ /μL	±0.5	1.5	1.5	1.5	1.67
MO	10 ³ /μL	±0.2	0.4	0.4	0.4	0.31
EO	10 ³ /μL	±0.2	0.1	0.1	0.1	0.15
BA	10 ³ /μL	±0.2	0.1	0.1	0.1	0.05
RBC	10 ⁶ /μL	±0.15	2.62	2.60	2.63	2.49
HGB	g/dL	±0.4	6.4	6.4	6.3	6.44
HCT	%	±2.0	21.6	20.8	20.5	22.2
MCV	fL	±7.0	82.5	80.0	77.9	89.3
MCH	pg	±2.4	24.4	24.6	24.0	25.9
MCHC	g/dL	±3.0	29.6	30.8	30.7	29.0
RDW-CV *	%	±4.0	14.9	16.2	15.2	15.2
PLT	10 ³ /μL	±30	50	57	51	50.1
MPV	fL	±3.0	4.1	5.2	4.4	6.9

Manual Registration	
Lot	248
Level	Low
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901690P