


**MEK-CAL** HEMATOLOGY CALIBRATOR 2.0mL

**ASSAY SHEET**

<b>LOT</b>	PLUS248
	2024-09-05
<b>OPEN</b>	7DAYS

MEK-9100,9200  
2D Barcode Entry



MEK-1301,1302,1303,1305  
Barcode Entry



		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 <sup>3</sup> /μL	9.2	9.2	9.4	9.3	9.1	9.47	9.19
RBC	10 <sup>6</sup> /μL	4.69	4.69	4.66	4.68	4.65	4.65	4.67
HGB	g/dL	13.5	13.6	13.5	13.4	13.3	13.65	13.65
HCT	%	40.1	40.1	41.8	40.5	39.5	45.7	46.0
MCV	fL	85.5	85.5	89.7	86.6	84.9	98.3	98.5
MCH	pg	28.8	29.0	29.0	28.6	28.6	29.4	29.2
MCHC	g/dL	33.7	33.9	32.3	33.1	33.7	29.9	29.7
PLT	10 <sup>3</sup> /μL	256	256	246	249	248	265.4	277.4
MPV	fL	7.6	9.3	7.6	8.0	7.6	9.8	9.9
RDW-CV*	%	15.6	15.6	14.6	15.8	14.9	15.8	15.7

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	167
FL	87
SD	108
FS THR.	100

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-901780R